

Dot-it-Spot-it Total Protein Assay

INSTRUCTIONS
Issued January 2025



Dot·it·
Spot·it



Dot-it-Spot-it Total Protein Kit.
 Store reagents refrigerated.
 Shipping at ambient temperature.

TABLE OF CONTENTS

	page
INTRODUCTION	3
REQUIRED ITEMS	3
WARRANTY	4
TECHNICAL ASSISTANCE	4
MANUFACTURE	4
PROCEDURE:	
Precautions	5
1. Preparations and sample buffers	5-6
2. Dispensing using a Template	7
3. Dispense and incubate	7
4. Detect and wash	8
5. Dry and mount	9
6. Scanning	10
7. Concentration estimation	10
NOTES:	
Note 1: Before starting	11-12
Note 2: Calibration curve, High samples	12
Note 3: Sample Buffer Compatibility Tables	13
Note 4: Scanning	16
Note 5: To get intensity values for each dot	17-19
Note 6: Concentration estimation	19-21

INTRODUCTION

- Dot-it-Spot-it measures the total protein concentration.
- One μL of sample mixture (application mixture) is applied on the Detection Sheet, which has positions for 20 samples.
- Sample proteins are binding to the Detection Sheet during a short incubation time.
- A black Detection Solution is flowed over the sheet and binds to the proteins, followed by a washing step.
- Finally the sheet is dried and mounted on a template, and the blackness intensity is quantified by the use of an image scanner.

REQUIRED ITEMS

Included in the Dot-it-Spot-it kit Store the reagents refrigerated. * The vial with Dilution Buffer needs to be warmed before use to dissolve the precipitates.	<ul style="list-style-type: none"> • Detection Sheet • Detection Solution • Washing Solution • Dilution Buffer* • Mounting Template
Equipment included in the Start-up package	<ul style="list-style-type: none"> • Light source • Dispensing Template with clips • Reservoirs with lids • Tape roller
Other equipment	<ul style="list-style-type: none"> • Pipettes • Hair dryer - or similar • Scissor
Detection equipment	<ul style="list-style-type: none"> • Image Scanner: Epson Perfection V39II (8-bit), or Perfection V600 Photo (16-bit)
Software recommendations	<ul style="list-style-type: none"> • ImageJ and Epson Scan 2
Files to download: dot-it-spot-it.com/method	<ul style="list-style-type: none"> • Scanning Template (xlsx) • Detection recommendations



Start-up package:

1. Reservoirs with lids
2. Dispensing Template with clips

WARRANTY

The information presented here is accurate to the best of our knowledge. It is the responsibility of the user to verify the suitability of the supplied materials and procedures for a particular purpose. In this respect further processing made by the user may affect the results, in which event Maplestone AB disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Maplestone AB and its authorized distributors, in such event, shall not be liable for damages indirect or consequential. The product is intended for research only.

TECHNICAL ASSISTANCE

If you experience any problems or difficulties with the products, please contact us.

Web: dot-it-spot-it.com
Email: info@dot-it-spot-it.com


MANUFACTURE

Maplestone AB, Högåsvägen 205, SE-74141 Knivsta, Sweden

PROCEDURE

Precautions Only handle the Detection Sheet by the mounted thicker pad. Do not touch the thin nitrocellulose sheet. It will directly capture proteins from your hands, traces on your gloves and any settled dust.

Even tiny drops of saliva can be seen on the sheet after a good laugh - it really is a protein sensitive technique! See Note 1.



Detection Sheet

Thick pad: You can touch and write on it.

Thin nitrocellulose sheet: A film of nitrocellulose on a polyester backing. Be careful, do not touch – it adsorbs proteins.

1. Preparations

Calibration standard

Bovine serum albumin (BSA) is often used as standard. Other proteins, or protein mixtures, with similar properties as the sample proteins might be preferred, see Note 2. It is important that the standard protein is in the same buffer as the samples. Prepare a dilution curve of the standard to each daily run.

Sample buffers

Dot-it-Spot-it assay can be used for a range of buffers containing both reducing substances and detergents. However, some additives in the sample buffer may restrict the percentage of sample you can use in the application mixture that you dispense to the Detection Sheet. Please refer to note 3 for a compatibility table featuring commonly used complex buffers and single buffer additives. In note 3, we refer to the compatibility for our Nano and Micro method.

Method	% sample	Application mixture	standard µg BSA/mL
Nano:	50	1 part of sample + 1 part DB	0.2-20
Micro:	5	1 part of sample + 50 part DB	2-200

Additionally, visit dot-it-spot-it.com/Results, Specificity and Compatibility, for updated information regarding sample buffer compatibility.

Application mixture

Take out the Dilution Buffer from the refrigerator and warm the vial until the precipitates are dissolved. Prepare the application mixture by mixing the sample/standard with the included Dilution Buffer (DB). Only 1 µL of application mixture is used for each dot replicate. The kit volume of Dilution Buffer is 15 mL, which is sufficient for using 50 µL per sample for dilution.

For a sample buffer that is fully compatible, such as TBS, you will achieve a more sensitive assay by increasing the amount of sample used in the application buffer. See the table below for details.

Sample in application mix	recommended low BSA standard µg/mL
50%	0,2
25%	0,4
10%	1,0
7%	1,4
5%	2,0
4%	2,5
2%	5
1%	10

2. Dispensing using a Template

1) Position the Detection Sheet on the Dispensing Template and use the clips to fasten it.



2) Place the light source under the Template and you will see the 20 circles through the thin nitrocellulose sheet.



3. Dispense and incubate

Dispense 1 μ L of your application mixture on the nitrocellulose sheet in the designed position shown by the template.

Incubate the Detection Sheet at room temperature for 30 min.

The required incubation time may vary due to the temperature, humidity and air circulation in the laboratory.



4. Detect and wash

This step requires the Detection and Washing Solutions. Invert the bottles before use. Fill up the two reservoirs, one with 1.5 mL of Detection Solution and the other with 1 mL of Washing Solution.

Take the Detection Sheet and press the tape and the thick pad together to make sure that liquid will be transported between nitrocellulose and the thick pad.



Place the nitrocellulose end of the Detection Sheet in the lid-covered reservoir, filled with the Detection Solution.



After 5 min, take up the Detection Sheet. Wipe off the surplus liquid at the end of the Detection Sheet. Move the Detection Sheet to the lid-covered reservoir, filled with the Washing Solution.



After 5 min, take up the Detection Sheet. Wipe off the surplus liquid at the end of the Detection Sheet.


The reservoirs are washed with water and dried before next use.




5. Dry and mount

Cut off the thicker pad using a scissor. Place the nitrocellulose sheet on the backside of the Dispensing Template.

Dry the nitrocellulose sheet using a hair dryer (or compressed air), to avoid migration of the colour in the dot.



Use the tape roller to mount the nitrocellulose sheet on the Mounting Template.



Make sure the nitrocellulose sheet is completely dry before scanning.

After scanning and calculating, dispose the nitrocellulose sheet (thin nitrocellulose film on a polyester backing) with dots of protein in accordance with your biohazard waste disposal protocol.


6. Scanning



Make sure that both scanner area and nitrocellulose sheet are free of dust.

Position the Mounting Template with the nitrocellulose sheet against the glass surface. We have found it worthwhile to check where we have the lowest imprecision over the scanning glass surface. Please see our web PDF file titled Detection Recommendations.

See Note 4, showing the required settings for Epson software



Scan the indicated area of the nitrocellulose sheet and save the tiff-file for further processing.

Follow Note 4 for quantifying the signal intensity of each dot using ImageJ and Excel.

7. Concentration estimation

Calculate the concentration for unknown samples by comparing their signal intensity values with the values for the calibration standard.

If plotting results by hand, a point-to-point fit is preferable to a linear fit.

Any software capable of four parameter logistic curve-fit, like the Rodbard function or similar, will work well for Dot-it-Spot-it.

A convenient curve-fitting can be obtained by using ImageJ – curve fitting, together with Excel, see Note 5.

NOTES

Note 1: Before starting

The Dot-it-Spot-it Total Protein Assay is an ultra sensitive test, using very low amounts of protein. When you are working with 0.1 nanogram protein you need to be aware of some pitfalls.

Make sure that you have clean tubes, clean pipette tips, clean gloves!

Plastic consumables, like tubes and tips, may be manufactured using slip agents and plasticisers that can leach into the sample. We are using Forensic DNA Grade Eppendorf tubes, and Eppendorf epT.I.P.S, with good results. The tubes and tips are handled with clean gloves, we order tips for 1 μ L already dispensed in trays. The immersion depth recommendation for tips are 1 mm for 1 μ L volume..

Do you generate additional proteins? If you place a finger on the nitrocellulose part of the Detection Sheet you may finally see the fingerprint. Saliva contains about 1 mg protein per mL, so be careful if talking or laughing. If you work with protein stock solutions of 1 mg/mL or higher, remember that the concentration is several thousand times higher than the concentration you measure, and your gloves might be contaminated. We work wearing a laboratory coat, disposable gloves, mob cap for covering the hair, and a medical mask.

Avoid protein adsorption and protease activity. For protein concentrations below 0.1 mg/mL there will be a substantial loss of protein to tube walls. Detergents, like Tween 20, are commonly added to reduce loss of protein. Presence of proteolytic enzymes can in some seconds degrade the proteins. Protease inhibitors, like EDTA, might be necessary to add.



Note 2:

Calibration curve

The need for a calibration curve (standard). Calibrate with a protein or protein mixture that corresponds to your samples.

The concentration of pure proteins can be checked by their absorbance at 280 nm, using 1 cm pathway in the photometer.

1 mg/mL BSA = 0.67 AU

1 mg/mL HSA = 0.53 AU


1 mg/mL IgG = 1.37 AU

Dilute your calibration protein in the same buffer as you have your samples in.

High samples

If you see very dark half-filled dots on the nitrocellulose sheet after testing, you might have too high protein concentration in your sample.

Below you can see the results for 1 and 10 μ g/mL BSA in the application mixture, which are OK to measure. Samples with 100-10.000 μ g/mL need to be diluted.



Note 3: Sample Buffer Compatibility Tables

In addition to proteins, sample buffers contain various essential additives for extracting, dissolving, and storing proteins from biological samples. These additives may affect the results of total protein assays.

The tables below indicate the recommended maximum percentage of sample in the application buffer when using a particular sample buffer. They also specify the recommended BSA concentration of the lowest standard point.

Table 1 displays the results for several complex buffers.

Table 2 presents some individual additives tested in TBS buffer with 1 mM EDTA. Additionally, except for the Detergent group, 0.03% Tween 20 was included. Both zero and BSA samples were analysed. The additive was considered satisfactory if the results with and without the additive were similar.

Please use the information in the tables as a reference only. It is the user's responsibility to confirm the compatibility of their specific sample buffer.

Table 1

	Application mixture	Recommended low standard $\mu\text{g BSA/mL}$
Complex buffers that can be directly applied		
<i>Method Nano-special</i>		
Laemmli-B buffer (type Biorad) mixed with PBS 33 mM TRIS-buffer, 5 mM phosphate, 75mM NaCl 13.2% glycerol, 1.1% SDS, 0.005% BFB	100% sample	0,15
Laemmli-S buffer (type Sigma) mixed with PBS 62,5 mM TRIS-buffer, 5 mM phosphate, 75 mM NaCl 10% glycerol, 2% SDS, 0.004% BFB	100% sample	0,15
OGP buffer mixed with PBS 10 mM TRIS-buffer, 16 mM phosphate, 150 mM NaCl 1% OGP, 1 mM EDTA	100% sample	0,15
Complex buffers that require further dilution with our Dilution Buffer		
<i>Method with 4-7 percentage of sample in the application mixture</i>		
RIPA buffer Pierce 25 mM TRIS-buffer, 150 mM NaCl 1% NP-40, 1% DOC, 0.1% SDS	4% sample	2,5
Cell Lysis buffer 25 mM TRIS-buffer, 150 mM NaCl 1% Triton X-100, 5% glycerol, 1 mM EDTA	7% sample	3
Buffer system that works well		
<i>Method Nano, Method Micro</i>		
TBS: 20 mM TRIS buffer pH 7.5 0.15 M NaCl, 0.02% NaN ₃	50% sample	0,2
	5% sample	2
PBS: 20 mM phosphate buffer pH 7.5 0.15 M NaCl, 0.02% NaN ₃	50% sample	0,2
	5% sample	2

Table 2

	Method Nano	Method Micro	unit
recommended low BSA standard	0.2 µg/mL	2 µg/mL	
Detergents			
Tween® 20	0,1	1,0	%
Triton™ X-100	0,1	1,0	%
SDS	4	20	%
Reducing agents			
TCEP	40	400	mM
DTT (fresh)	20	200	mM
BME	5	50	%
Sugars			
Lactose	40	400	mM
GlcNAc	20	200	mM
Protease inhibitors			
cOmplete™, EDTA free, Roche <i>w.s.= working solution</i>	0.5 x	1 x	w.s.
EDTA	1	1	mM
PMSF	1	1	mM
Buffer additives			
Glycerol	5	50	%
EDTA	20	200	mM
urea	6	10	M
thiourea	2	2	M

Abbreviations

- BFB: Bromophenol blue
- BSA: Bovine serum albumin
- BME: 2-Mercaptoethanol, beta-Mercaptoethanol
- DOC: Deoxycholic acid
- DTT: Dithiothreitol
- EDTA: Ethylenediaminetetraacetic acid
- GlcNAc: N-Acetylglucosamine
- PMSF: Phenylmethylsulfonyl fluoride
- SDS: Sodium dodecyl sulphate
- TCEP: Tris(2-carboxyethyl)phosphine


Note 4: Scanning

Instructions for using the software Epson Scan 2 or SilverFast.

Ensure correct settings in the scan software:

- Document type: Photo,
- Output image: tiff-file.
- Bit depth: Grayscale16-bit or 8-bit.
- Pixel resolution: 600 dpi.
- Adjustments: It is important that the image is not adjusted by the software.
- Configuration/Colour: No colour correction.
- Gamma gradation: Epson scan 1.8, SilverFast 1.0.

If you like to have an overview image , save a pdf file of the complete Mounting Template..



Note 5: To get intensity values for each dot

Download the free software ImageJ: <https://imagej.net/ij/download.html>




Download the Excel file named Scanning Template: Set-it-spot-it.com/method

Set the conditions for ImageJ measurement: *Analyze / Set measurement.*

For 16-bit value, use measurement of Median and Mean gray value.

For 8-bit values, use Decimal places= 3, and Mean gray value.



Open your tiff image of the nitrocellulose strip in Image J.



Use the oval tool and measure the dot - use "Analyze/ Measure" or use the rapid command "M". An ImageJ results table will collect the values.

We recommend to measure the whole dot area, but not outside this area, as shown in the image. However, if you see tiny black spots inside the dot you can limit the measurement area to avoid the small contamination spots.



Reflected light values. The values for the 20 dots are collected in an ImageJ results table. Reflected light: a dark dot gives low value, whilst a white dot gives a high value.

	Mean	Median
1	63102	63157
2	62091	62312
3	60653	60653
4	55759	55764
5	46209	46277
6	34102	33861
7	62950	63032
8	62295	62347
9	60573	60542
10	55364	55366
11	46011	46061
12	34188	34138
13	62834	62877
14	62214	62262
15	60393	60451
16	55652	55665
17	46002	46080
18	33613	33637

	Mean
1	235.856
2	239.266
3	236.588
4	228.791
5	208.538
6	170.611
7	117.817
8	236.480
9	234.634
10	227.945
11	206.854
12	168.307
13	114.377
14	233.236
15	234.401
16	228.945
17	210.578
18	170.549
19	114.937
20	230.770

The 16-bit median light intensity gives values in the range 0–65535. For 8 bit values, use the mean value with three decimal places for the light intensity, which lies in the range 0–251,000.

Blackness per pixel values. Copy the Image J results table into worksheet tab "1. Paste from ImageJ" in the Excel file Scanning Template. The median blackness per pixel (bpp) will be calculated as RESULTS. Blackness per pixel: a dark dot will have a higher bpp value than a light dot.

Dot-It-Spot-It Protein Assay		Blackness per pixel (bpp), 16-bit greyscale= 0, 1-65.535 levels					Reflected light	
v. 210202		Blackness per pixel = bpp					PASTE FROM IMAGEJ	
Detection SHEET	Sample Position	RESULTS Median bpp	bpp= 65.535-Reflected light			Mean	Median	
			Median bpp	Mean bpp	Mean/Median	Mean	Median	
1	1	2 378,0	1	2 378,0	2 433,0	1,02	63102	63157
1	2	3 223,0	2	3 223,0	3 444,0	1,07	62091	62312
1	3	4 882,0	3	4 882,0	4 882,0	1,00	60653	60653
1	4	9 771,0	4	9 771,0	9 776,0	1,00	55759	55764
1	5	19 258,0	5	19 258,0	19 326,0	1,00	46209	46277
1	6	31 674,0	6	31 674,0	31 433,0	0,99	34102	33861
1	7	2 503,0	7	2 503,0	2 585,0	1,03	62950	63032
1	8	3 188,0	8	3 188,0	3 240,0	1,02	62295	62347
1	9	4 993,0	9	4 993,0	4 962,0	0,99	60573	60542
1	10	10 169,0	10	10 169,0	10 171,0	1,00	55364	55366
1	11	19 474,0	11	19 474,0	19 524,0	1,00	46011	46061
1	12	31 397,0	12	31 397,0	31 347,0	1,00	34188	34138
1	13	2 658,0	13	2 658,0	2 701,0	1,02	62834	62877
1	14	3 273,0	14	3 273,0	3 321,0	1,01	62214	62262
1	15	5 084,0	15	5 084,0	5 142,0	1,01	60393	60451
1	16	9 870,0	16	9 870,0	9 883,0	1,00	55652	55665
1	17	19 455,0	17	19 455,0	19 533,0	1,00	46002	46080
1	18	31 898,0	18	31 898,0	31 922,0	1,00	33613	33637

The 16-bit median blackness per pixel (bpp) gives test values in the range 2,000–50,000 bpp.

Dot-It-Spot-It Protein Assay calculation template for 8 bit greyscale
v. 241103
Blackness per pixel = bpp
bpp = 255-Reflected light

RESULTS		Reflected light	
Detection SHEET	Sample Position	Mean	
		PASTE FROM IMAGEJ	
		Mean	
1	1	24,079	230,921
1	2	15,747	239,253
1	3	18,494	236,506
1	4	26,409	228,591
1	5	46,766	208,234
1	6	84,803	170,197
1	7	137,251	117,749
1	8	18,675	236,325
1	9	20,348	234,652
1	10	27,041	227,959
1	11	48,101	206,899
1	12	86,785	168,215
1	13	140,039	114,961
1	14	32,223	222,777
1	15	20,540	234,46
1	16	25,938	229,062
1	17	44,036	210,964
1	18	84,263	170,737
1	19	139,444	115,556
1	20	24,181	230,819

The mean 8-bit blackness per pixel (bpp) gives test values in the range 12–150 bpp.

Note 6: Concentration estimation for unknown samples

A convenient curve-fitting and concentration estimation can be obtained by combining ImageJ – curve fitting, and the Excel file Scanning Template, worksheet tab “3. Concentration”.

Calculate the mean bpp value for your standard replicates. Worksheet tab “2 Standard” might be useful, remember that ImageJ use dots and not comma. Curve fitting ImageJ: see Analyze – Tools – Curve Fitting, and choose Rodbard.

The first table row is for concentration values for your standard points, and the second for the obtained bpp values. Fill in your values and push Fit.



Rodbard	
0.00	2378
0.1	3223
0.3	4882
1.0	9771
3.0	19258
10.0	31674

Rodbard	
0	16.9255
0.1	18.3105
0.3	24.5915
1.0	48.41
3.0	88.292
10	144.339

Test values for µg BSA/mL standard in 16-bit (left Fig.) and 8-bit bpp.

Parameters from ImageJ

After fitting, you can see how the individual standard concentrations forms a standard curve. The parameters a-d will be used in Template of Calculation.



Standard curve for µg BSA/mL (x-axis) with 16-bit (upper Fig.) and 8-bit bpp signals.

Paste parameters a-d from ImageJ into the Excel file.

In Scanning Template, worksheet tab "3. Concentration", you can find the signal values from tab "1. Paste from ImageJ". Fill in the values for parameter a-d, and the concentration for the unknown samples will be calculated.

Calculation of concentration using curve fit Rodbard

Fill in the parameters obtained from ImageJ for the dilution curve va

Rodbard curve fit	
a	
b	
c	
d	

Fill in values obtained from ImageJ

Detection		SIGNAL	Conc.	
Sheet	Sample position	[dbpp]	[ug/mL]	
1	1	2378	#####	Standard
1	2	3223	#####	Standard

Calculation of concentration using curve fit Rodbard

Fill in the parameters obtained from ImageJ for the dilution curve va

Rodbard curve fit	
a	2404
b	1,04
c	4,35
d	43991

Fill in values obtained from ImageJ

Detection		SIGNAL	Conc.	
Sheet	Sample position	[dbpp]	[ug/mL]	
1	1	2378	#OGILTGT!	Standard
1	2	3223	0,10	Standard
1	3	4882	0,31	Standard
1	4	9771	0,99	Standard
1	5	19258	3,01	Standard
1	6	31674	10,00	Standard
1	7	2503	0,01	Samples
1	8	3188	0,10	Samples
1	9	4993	0,32	Samples
1	10	10169	1,06	Samples
1	11	19474	3,07	Samples
1	12	31397	9,70	Samples
1	13	2658	0,03	Samples
1	14	3273	0,11	Samples
1	15	5084	0,33	Samples
1	16	9870	1,01	Samples
1	17	19455	3,07	Samples